

# Skytron Elite 3600B Premium Series Replacement Pad Set



*Stock Photo - Options May Vary*

## FACTORY NEW

The Skytron Elite 3600B Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of medium-firm high-density foam combined with 1 inch of Visco memory foam, providing optimal support and pressure distribution. The conductive 4-way stretch Fusion III soft vinyl cover ensures durability and ease of cleaning in the operating room environment.

## FEATURES

- Designed for Longer-Duration Procedures
- Hook & Loop Kit comes with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)