

Skytron 3100 Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 3100 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2 inches of medium-firm high-density foam construction covered with electrically conductive vinyl for reliable performance in the operating room.

This replacement pad set includes a hook and loop kit for secure attachment to the surgical table. The waterfall edge construction and fluid-resistant cover help maintain the pad's integrity during surgical procedures.

The pad is X-ray permeable, anti-microbial, and latex-free, making it suitable for a wide range of surgical applications and patient populations.

FEATURES

- Designed for Short-Duration Procedures
- Hook & Loop Kit comes included with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)