

# Skytron Elite 3500 Standard Series Replacement Pad Set



*Stock Photo - Options May Vary*

## FACTORY NEW

The Skytron Elite 3500 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. This pad set is constructed with medium-firm, high-density foam that is 2 inches thick and covered with electrically conductive vinyl for safety and durability. The set is compatible with both the Skytron Elite 3500 and Elite 3500B surgical tables.

## FEATURES

- Designed for Short-Duration Procedures
- Hook & Loop Kit comes with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)