


Skytron Elite 3500B Standard Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 3500B Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2 inches of medium-firm high-density foam construction covered with electrically conductive vinyl for safety and durability.

This replacement pad set is compatible with both the Skytron Elite 3500B and Elite 3500 models. The pad includes a hook and loop kit for secure attachment to the surgical table.

For bariatric patients, longer procedures, or enhanced comfort requirements, premium pad options are available as alternatives to this standard series.

FEATURES

- Designed for Short-Duration Procedures
- Hook & Loop Kit comes with Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)