

Skytron Elite 3500B Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 3500B Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of medium-firm high-density foam layered with 1 inch of visco memory foam, providing optimal support and pressure distribution. The set is compatible with both the Elite 3500B and Elite 3500 models.

The pad surface is covered with conductive 4-way stretch Fusion III soft vinyl, designed to withstand the demanding surgical environment. The waterfall edge construction ensures secure positioning and durability throughout the pad's service life.

This replacement pad set includes a hook and loop kit for easy installation and attachment to the surgical table base.

FEATURES

- Designed for Longer-Duration Procedures
- Hook & Loop Kit comes with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)