

# Skytron 3501 Standard Series Replacement Pad Set



*Stock Photo - Options May Vary*

## FACTORY NEW

The Skytron 3501 Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. The pad features 2" thick medium-firm high-density foam construction covered with electrically conductive vinyl for reliable performance during surgical applications.

This replacement pad set is compatible with both the Skytron 3501 and 3501B surgical tables. The pad includes a Velcro kit for secure attachment and is constructed with waterfall edge design for optimal patient positioning and comfort during procedures.

## FEATURES

- Designed for Short-Duration Procedures
- Velcro Kit Comes with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)