

Skytron Elite 3501 Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 3501 Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of medium-firm high-density foam combined with 1 inch of Visco memory foam, providing optimal support and pressure distribution. The set is compatible with both the Skytron Elite 3501 and Elite 3501B surgical tables.

The pad surface is covered with conductive 4-way stretch Fusion III soft vinyl that resists fluid and stains while maintaining X-ray permeability. The waterfall edge construction ensures secure positioning and durability. Each replacement pad set includes a hook and loop attachment kit for easy installation.

The pad is constructed from latex-free materials with anti-microbial properties to support infection control protocols in the operating room.

FEATURES

- Designed for Longer-Duration Procedures
- Hook & Loop Kit comes with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)