

# Skytron 3501B Standard Series Replacement Pad Set



Stock Photo - Options May Vary

## FACTORY NEW

The Skytron 3501B Standard Series Replacement Pad Set is designed for low-risk patients and short-duration procedures. This pad set is constructed with medium-firm high-density foam and covered with electrically conductive vinyl, providing a durable and reliable surface for surgical applications. The set is compatible with both the Skytron 3501B and 3501 models.

The pad features a 2-inch thickness and waterfall edge construction for optimal patient positioning and comfort. It is fluid and stain resistant, making it easy to maintain and clean between procedures. The X-ray permeable design allows for imaging without pad removal.

The complete replacement set includes a velcro kit for secure attachment to the surgical table. The pad is latex-free and treated with anti-microbial protection to support infection control protocols in the operating room.

## FEATURES

- Designed for Short-Duration Procedures
- Velcro Kit comes Complete with Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

## SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)