

Skytron Elite 3501B Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron Elite 3501B Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features a 2-inch overall thickness constructed with medium-firm high-density foam and 1-inch Visco memory foam, providing optimal support and pressure relief.

The pad is covered with conductive 4-way stretch Fusion III Soft Vinyl, designed to withstand the demanding environment of the operating room. This replacement pad set is compatible with both the Skytron Elite 3501B and Elite 3501 surgical tables.

FEATURES

- Designed for Longer-Duration Procedures
- Hook & Loop Kit comes with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4-Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)