

Skytron 6500HD Premium Series Replacement Pad Set



Stock Photo - Options May Vary

FACTORY NEW

The Skytron 6500HD Premium Series Replacement Pad Set is engineered for patient comfort during extended surgical procedures. The pad features 2 inches of medium-firm high-density foam with 1 inch of Visco memory foam, covered in conductive 4-way stretch Fusion III soft vinyl for durability and ease of cleaning.

This replacement pad set is compatible with Skytron models 5001, 6001, 6001K, 6001KR, 6001KRB, 6500, 6500B, and 6500HD. The complete set includes a Velcro kit for secure attachment to the surgical table.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)