

Hill Rom TranStar - Replacement Stretcher Pad - 26" x 75" x 4"

FACTORY NEW

The Hill Rom TranStar replacement stretcher pad measures 26" wide by 75" long with a 4" thickness. It features three layers of high-density, high-resilient medical-grade urethane foam that provides outstanding pressure redistribution, postural support, and comfort for patients.

The pad is constructed with universal corners featuring two round and two angled cuts (foot-radius and head-angle) for versatile compatibility. It includes a fluid protection zipper guard to prevent contamination and meets fire safety standards with CA 117 rating.

The pad offers multiple protective properties including anti-microbial, anti-static, stain-resistant, abrasion-resistant, and non-allergenic characteristics. Loop Velcro is installed down the center bottom unless otherwise specified.



Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)