

Standard Stretcher Pad - 26" x 75" x 4"



Stock Photo - Options May Vary

FACTORY NEW

The standard stretcher pad measures 26" x 75" x 4" and is constructed with quality medium-density foam and a latex-free, anti-microbial, fluid-resistant cover. This combination ensures patient comfort while meeting the requirements of medical professionals.

These durable pads are recommended for minimally invasive procedures or lower-risk patients. The design accommodates many styles and sizes of stretchers, making them a versatile choice for various clinical settings.

The pad features easy-to-clean construction and is built with high-quality materials designed for reliable performance in medical environments.

FEATURES

- 26" x 75" x 4" Thick
- Easy to clean
- Antimicrobial
- Latex-free and fluid-resistant
- Pressure-management padding available
- High-quality material and construction

SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)