

Stryker 1210 - Replacement Stretcher Pad - 26" x 76" x 4"

FACTORY NEW

The Stryker 1210 replacement stretcher pad measures 26" wide by 76" long with a 4" thickness, designed as a universal fit for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. The pad features universal corners with two round and two angled cuts, a fluid protection zipper guard, and is treated to be anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic.



Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)