

Stryker 1501 - Replacement Stretcher Pad - 29" x 76" x 4"

FACTORY NEW

The Stryker 1501 replacement stretcher pad measures 29" wide by 76" long with a 4" thickness, designed as a universal fit for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. The pad features universal corner design with two round and two angled cuts, fluid protection with zipper guard, and is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic for reliable performance in clinical environments.



Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 29"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)