

# Stryker 1509 - Replacement Stretcher Pad - 29" x 76" x 4"

## FACTORY NEW

The Stryker 1509 replacement stretcher pad measures 29" wide by 76" long with a 4" thickness, designed as a universal fit for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. The pad features universal corners with two round and two angled cuts, a fluid protection zipper guard, and is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic for durability and hygiene in clinical environments.



*Stock Photo - Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 29"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)