

Stryker 1731 - Replacement Stretcher Pad - 29" x 76" x 4"

FACTORY NEW

The Stryker 1731 is a replacement stretcher pad constructed from three layers of high-density, high-resilient medical-grade urethane foam. This construction provides outstanding pressure redistribution, postural support, and comfort for patients during transport and care. The pad measures 29 inches wide by 76 inches long with a 4-inch thickness, making it suitable for standard stretcher applications.



The pad features universal corners with two round and two angled cuts to accommodate various stretcher designs. It includes a fluid protection zipper guard to prevent contamination and is equipped with loop Velcro along the center bottom for secure attachment unless otherwise specified.

The pad is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic, making it suitable for healthcare environments. It meets fire safety standards with CA 117 rating and weighs 25 pounds.

Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 29"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)