

Steris-Birkova Stretcher Patient 3" All Velcro Safety Straps: PR-9200



Stock Photo - Options May Vary

FACTORY NEW

The Steris-Birkova Stretcher Patient 3" Velcro Safety Straps PR-9200 deliver reliable support for patients during transport. The straps work seamlessly with standard medical stretchers to secure patients quickly, improving workflow during emergencies and reducing delays. The strap design provides patients with a greater sense of stability, helping reduce stress during transfers.

The PR-9200 measures 3 inches in width and uses durable Dura-Cloth material for consistent performance in demanding environments. Velcro closures grip firmly and allow fast adjustments without slowing clinical routines. The straps clean easily, enabling staff to maintain hygiene standards with minimal effort.

Each PR-9200 set includes three two-piece straps that provide flexible placement across the stretcher. This configuration improves positioning control and helps medical teams adapt quickly to different patient needs.

FEATURES

- Made with Dura-Cloth
- All Velcro Straps easily mount on the stretcher
- Velcro Closure to secure/release quickly
- Helps maintain patient safety while on the stretcher
- Provides patients with a sense of security

SPECIFICATIONS

- Construction: Dura-Cloth
- Mounting: Metal Rings
- Closure: Velcro
- Straps are two individual straps and three pairs per set
- Dimensions per pair (2 Straps): 1 - 27.5" L x 3" W and 1 - 51" L x 3" W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)