

# Standard Stretcher Pad - 24" x 75" x 4"

## FACTORY NEW

The Standard Stretcher Pad measures 24" x 75" x 4" and is constructed with quality medium-density foam and a latex-free, anti-microbial, fluid-resistant cover. The pad is designed to ensure patient comfort while meeting the requirements of medical professionals.

These durable stretcher pads are recommended for minimally invasive procedures or lower-risk patients. The universal design accommodates many styles and sizes of stretchers commonly used in North American healthcare facilities.

The pad features easy-to-clean construction and comes with universal corners (2 round and 2 angled) for compatibility with various stretcher models. Loop Velcro is positioned down the center bottom unless otherwise specified.



*Stock Photo - Options May Vary*

## FEATURES

- 24" x 75" x 4" Thick
- Easy to clean
- Antimicrobial
- Latex-free and fluid-resistant
- High-quality material and construction
- Universal: fit most North American Stretchers
- Universal Corners: 2 Round and 2 Angled (Foot-Radius Cut, Head-angle cut)
- Custom Sizing is Available

## SPECIFICATIONS

- Width: 24"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)