

Stryker 1105 Prime - Replacement Stretcher Pad - 26" x 75" x 4"



Stock Photo - Options May Vary

FACTORY NEW

The Stryker 1105 Prime Replacement Stretcher Pad measures 26" wide by 75" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad features three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. Universal corner design includes two round corners and two angled corners (foot-radius cut and head-angle cut) to accommodate various stretcher models.

The pad incorporates a fluid protection zipper guard to prevent contamination from bodily fluids and other liquids. The cover material is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic, making it suitable for high-use clinical environments. The pad is fire-rated to CA 117 standards for safety compliance.

Loop Velcro is installed down the center bottom unless otherwise specified. The pad weighs 25 lbs and is available in multiple configuration options to meet specific facility requirements.

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)