

Stryker 720 – Replacement Stretcher Pad – 24" x 76" x 4"

FACTORY NEW

The Stryker 720 replacement stretcher pad is constructed from three layers of high-density, high-resilient medical-grade urethane foam. This construction provides outstanding pressure redistribution, postural support, and comfort for patients during transport and care. The pad measures 24 inches wide by 76 inches long with a 4-inch thickness, making it suitable for standard stretcher applications.



The pad features universal corners with two round corners at the head and two angled corners at the foot for secure fit and positioning. A fluid protection zipper guard protects against contamination, while the surface is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic. The pad is fire-rated to CA 117 standards for safety compliance.

Loop Velcro is installed down the center bottom unless otherwise specified. The durable construction ensures long-term performance in demanding medical environments.

Stock Photo – Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 24"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)