

Stryker 721 – Replacement Stretcher Pad – 24" x 76" x 4"

FACTORY NEW

The Stryker 721 replacement stretcher pad measures 24" wide by 76" long with a 4" thickness, designed as a universal fit for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. The pad features universal corners with two round and two angled cuts, loop Velcro down the center bottom, and a fluid protection zipper guard to prevent contamination.



Stock Photo – Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 24"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)