

Standard Stretcher Pad - 27" x 75" x 4"

FACTORY NEW

Standard stretcher pads are constructed with quality medium-density foam and a latex-free, anti-microbial, fluid-resistant cover designed to ensure patient comfort while meeting medical professionals' needs. These durable pads are recommended for minimally invasive procedures or lower-risk patients and can accommodate many styles and sizes of stretchers. The pads feature universal corners with 2 round and 2 angled cuts to fit most North American stretcher models.



Stock Photo - Options May Vary

FEATURES

- Dimensions: 27" x 75" x 4" Thick
- Easy to clean
- Pressure-management padding available
- Universal Corners: 2 Round and 2 Angled (Foot-Radius Cut, Head-angle cut)
- Universal: Fits most North American Stretchers
- Custom Sizing Available

SPECIFICATIONS

- Width: 27"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)