

# Standard Stretcher Pad - 27" x 77" x 4"

## FACTORY NEW

Standard stretcher pad measuring 27" wide by 77" long with 4" thickness. Constructed with quality medium-density foam and a latex-free, anti-microbial, fluid-resistant cover designed for patient comfort and durability.

Recommended for minimally invasive procedures and lower-risk patients. The pad features universal corners with 2 round and 2 angled cuts to accommodate many styles and sizes of stretchers commonly used in North American healthcare facilities.

The pad is easy to clean and maintain, with loop Velcro positioned down the center bottom for secure attachment. High-quality construction ensures reliable performance in clinical settings.



*Stock Photo - Options May Vary*

## FEATURES

- 27" x 77" x 4" Thick
- Universal Corners: 2 Round and 2 Angled (Foot-Radius Cut, Head-angle cut)
- Universal: fit most North American Stretchers
- Easy to clean
- Pressure-management padding available
- High-quality material and construction
- Electrically Conductive
- Custom Sizing Available

## SPECIFICATIONS

- Width: 27"
- Length: 77"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)