

# Standard Stretcher Pad - 30" x 76" x 4"

## FACTORY NEW

Standard stretcher pad measuring 30" wide by 76" long with 4" thickness, constructed from quality medium-density foam. The pad features a latex-free, anti-microbial, and fluid-resistant cover designed to ensure patient comfort while meeting medical professionals' requirements.

Recommended for minimally invasive procedures and lower-risk patients, this stretcher pad accommodates many styles and sizes of stretchers. Universal corners feature 2 round and 2 angled cuts to fit various stretcher configurations.

The pad includes loop Velcro down the center bottom for secure positioning. Custom sizing options are available to meet specific facility needs.



*Stock Photo - Options May Vary*

## FEATURES

- X-Ray Permeable
- Universal: fits most North American Stretchers
- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Resistant / Flame Retardant
- Electrically Conductive
- Custom Sizing Available

## SPECIFICATIONS

- Width: 30"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)