

Stryker 1115 Prime - Replacement Stretcher Pad - 30" x 75" x 4"

FACTORY NEW

The Stryker 1115 Prime Replacement Stretcher Pad measures 30" x 75" x 4" and is constructed with three layers of high-density, high-resilient medical-grade urethane foam. This multi-layer construction provides outstanding pressure redistribution, postural support, and patient comfort.



The pad features universal corners with two round and two angled cuts (foot-radius and head-angle) to accommodate various stretcher models. A fluid protection zipper guard protects against fluid contamination, while the surface is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic.

The pad is fire-rated to CA 117 standards and comes with loop Velcro down the center bottom unless otherwise specified. This durable stretcher pad is designed for reliable performance in clinical settings.

Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 30"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)