

# Stryker 1115 Prime - Replacement Stretcher Pad - 26" x 75" x 4"

## FACTORY NEW

The Stryker 1115 Prime Replacement Stretcher Pad measures 26" wide by 75" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad features three layers of high-density, high-resilient medical-grade urethane foam that provides outstanding pressure redistribution, postural support, and patient comfort. Universal corner design includes two round corners and two angled corners (foot-radius cut and head-angle cut) to accommodate various stretcher configurations.

The pad incorporates a fluid protection zipper guard to prevent fluid contamination and maintain a hygienic patient surface. It meets fire safety standards with CA 117 rating and offers comprehensive protective properties including anti-microbial, anti-static, stain-resistant, abrasion-resistant, and non-allergenic characteristics. Loop Velcro is positioned down the center bottom for secure attachment unless otherwise specified.



*Stock Photo - Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)