

Stryker 1105 Prime - Replacement Stretcher Pad - 30" x 75" x 4"

FACTORY NEW

The Stryker 1105 Prime Replacement Stretcher Pad measures 30" x 75" x 4" and is constructed with three layers of high-density, high-resilient medical-grade urethane foam. This multi-layer construction provides outstanding pressure redistribution, postural support, and patient comfort during transport and care.



The pad features universal corners with two round and two angled cuts designed for proper fit and positioning. A fluid protection zipper guard protects against contamination, while the surface is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic.

The pad meets fire safety standards (CA 117) and comes with loop Velcro down the center bottom for secure attachment to stretcher frames. The durable construction ensures reliable performance in demanding medical environments.

Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 30"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)