

# Stryker 737 – Replacement Stretcher Pad – 26" x 75" x 4"

## FACTORY NEW

The Stryker 737 replacement stretcher pad measures 26" wide by 75" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad features three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. Universal corner design with rounded corners at the foot and angled corners at the head accommodates various stretcher models.

The pad incorporates a fluid protection zipper guard to prevent contamination from bodily fluids and other liquids. The cover material is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic, making it suitable for frequent cleaning and high-use clinical environments.

The pad is fire-rated to meet CA 117 standards and comes with loop Velcro positioned down the center bottom for secure attachment to stretcher frames. The durable construction supports extended use in hospitals, ambulances, and other medical transport applications.



*Stock Photo – Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)