

Stryker 738 – Replacement Stretcher Pad – 30" x 75" x 4"

FACTORY NEW

The Stryker 738 replacement stretcher pad measures 30" wide by 75" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad features three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. Universal corners with rounded and angled cuts accommodate various stretcher configurations, while loop Velcro along the center bottom ensures secure positioning.



Stock Photo – Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 30"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)