

Stryker 735 – Replacement Stretcher Pad – 26" x 75" x 4"

FACTORY NEW

The Stryker 735 replacement stretcher pad is constructed with three layers of high-density, high-resilient medical-grade urethane foam designed to provide outstanding pressure redistribution, postural support, and comfort. The pad measures 26 inches wide by 75 inches long with a 4-inch thickness, making it suitable for standard stretcher applications.



This durable pad features universal corners with two round and two angled cuts to accommodate various stretcher designs. It includes a fluid protection zipper guard to prevent contamination and is treated to be anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic. The pad is fire rated to CA 117 standards for safety compliance.

Stock Photo – Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)