

# Stryker Atlas 660 – Replacement Stretcher Pad – 31" x 76" x 4"

## FACTORY NEW

The Stryker Atlas 660 replacement stretcher pad measures 31" wide by 76" long with a 4" thickness, designed as a direct replacement for standard Atlas stretcher systems. The pad features three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort during transport and care.



*Stock Photo – Options May Vary*

The pad incorporates universal corners with two round and two angled cuts (foot-radius and head-angle) to fit various stretcher configurations. A fluid protection zipper guard protects against contamination, while the surface is anti-microbial, anti-static, stain-resistant, abrasion-resistant, and non-allergenic for enhanced durability and hygiene.

The pad meets fire safety standards (CA 117 rated) and comes with loop Velcro along the center bottom for secure attachment unless otherwise specified. The construction ensures long-lasting performance in demanding medical environments.

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 31"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)