

Hill Rom TranStar - Replacement Stretcher Pad - 30" x 75" x 4"

FACTORY NEW

The Hill Rom TranStar replacement stretcher pad measures 30" x 75" x 4" and is constructed with three layers of high-density, high-resilient medical-grade urethane foam. This multi-layer construction provides outstanding pressure redistribution, postural support, and patient comfort.

The pad features universal corners with two round and two angled cuts (foot-radius and head-angle) to accommodate various stretcher models. A fluid protection zipper guard ensures protection against fluid contamination during use and cleaning.

The pad is anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic. It meets fire safety standards (CA 117) and comes with loop Velcro down the center bottom as standard unless otherwise specified.



Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 30"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)