

Hausted 462 Series - Replacement Stretcher Pad - 30" x 77" x 4"



Stock Photo - Options May Vary

FACTORY NEW

The Hausted 462 Series replacement stretcher pad measures 30" wide by 77" long with a 4" thickness, designed as a direct replacement for the Hausted 462 Horizon Series. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort.

The pad features universal corners with two round and two angled cuts (foot-radius and head-angle) to accommodate various stretcher configurations. It includes a fluid protection zipper guard to prevent contamination and is treated to be anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic.

The pad meets fire safety standards (CA 117 rated) and comes with loop Velcro positioned down the center bottom unless otherwise specified. The durable construction is engineered to withstand the demands of healthcare environments while maintaining patient comfort and safety.

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 30"
- Length: 77"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)