

Stryker 1211 - Replacement Stretcher Pad - 26" x 76" x 4"

FACTORY NEW

The Stryker 1211 replacement stretcher pad is constructed from three layers of high-density, high-resilient medical-grade urethane foam. This construction provides outstanding pressure redistribution, postural support, and comfort for patients during transport and care.

The pad measures 26 inches wide by 76 inches long with a 4-inch thickness, making it compatible with standard stretcher systems. Universal corners feature two rounded corners at the head and two angled corners at the foot for secure fit and positioning.

The pad is engineered with anti-microbial, anti-static, stain-resistant, and abrasion-resistant properties. It includes a fluid protection zipper guard to prevent contamination and is fire-rated to CA 117 standards.



Stock Photo - Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)