

# Stryker 1009 - Replacement Stretcher Pad - 26" x 76" x 4"

## FACTORY NEW

The Stryker 1009 is a replacement stretcher pad constructed from three layers of high-density, high-resilient medical-grade urethane foam. This construction provides outstanding pressure redistribution, postural support, and comfort for patients during transport and care.

The pad measures 26 inches wide by 76 inches long with a 4-inch thickness, making it suitable for standard stretcher applications. Universal corners feature two round and two angled cuts to accommodate different stretcher designs.

The pad is engineered with anti-microbial, anti-static, stain-resistant, and abrasion-resistant properties. It includes a fluid protection zipper guard to prevent contamination and meets fire safety standards.



*Stock Photo - Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 26"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)