

Stryker 721 – Replacement Stretcher Pad – 26" x 76" x 4"

FACTORY NEW

The Stryker 721 replacement stretcher pad measures 26" wide by 76" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilience medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. Its durable construction resists fluid contamination, staining, abrasion, and microbial growth while meeting fire safety standards.



Stock Photo – Options May Vary

FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

SPECIFICATIONS

- Width: 26"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)