

# Stryker 720 – Replacement Stretcher Pad – 26" x 76" x 4"

## FACTORY NEW

The Stryker 720 replacement stretcher pad measures 26" wide by 76" long with a 4" thickness, designed as a direct replacement for standard stretcher systems. The pad is constructed from three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort. The pad features universal corners with two round and two angled cuts, loop Velcro down the center bottom, and is available in various configurations to meet specific facility requirements.



*Stock Photo – Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 26"
- Length: 76"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)