

# Stryker 1005/SM104 - Replacement Stretcher Pad - 26" x 75" x 4"

## FACTORY NEW

The Stryker 1005/SM104 (M-Series) replacement stretcher pad measures 26" wide by 75" long with a 4" thickness, designed to fit standard stretcher frames. The pad is constructed from three layers of high-density, high-resilient medical-grade urethane foam that provides excellent pressure redistribution, postural support, and patient comfort.

The pad features universal corners with two round and two angled cuts (foot-radius and head-angle) for secure fit. It includes a fluid protection zipper guard to prevent contamination and is treated to be anti-microbial, anti-static, stain resistant, abrasion resistant, and non-allergenic.

The pad is fire-rated to CA 117 standards and comes with loop Velcro down the center bottom unless otherwise specified. The durable construction ensures long-term performance in clinical settings.



*Stock Photo - Options May Vary*

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117

## SPECIFICATIONS

- Width: 26"
- Length: 75"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)