

Blickman (Snyder) Solution Stand | Model 7808SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7808SS - (Snyder) Solution Stand. Featuring a double basin setup, providing greater capacity and functionality for handling multiple solutions or tasks simultaneously.

The basins are supported by a robust frame, which is further reinforced by a bottom shelf, offering additional storage space for equipment or supplies. Constructed entirely from stainless steel with welded joints, the stand ensures superior strength and stability, making it ideal for demanding environments where reliability is crucial.

The stainless steel construction also ensures easy maintenance and long-lasting performance.

FEATURES

- Stainless Steel Construction
- Double Basin with Supporting Frame
- Supporting Frame Connected by a Bottom Shelf
- Enhanced Strength and stability
- Large Basin Capacity

SPECIFICATIONS

- Construction: Seamless, Welded, Stainless Steel
- Stainless Steel H-Bracket
- 2" Swivel Casters
- Dimensions: 29"L x 15"W x 33"H
- Approx. Shipping Weight: 33 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)