

Blickman (Snyder) Solution Stand | Model 7808SS-HB



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7808SS-HB - (Snyder) Solution Stand: Featuring a double basin design, providing the capability to handle multiple tasks simultaneously with ease. The basins are supported by a robust frame with an H-brace, providing superior stability and strength for heavy usage.

Each stainless steel basin has a generous 8-1/2 quart capacity, ideal for larger volumes of solutions or materials. Despite their large capacity, the basins remain unusually quiet when dropped into the circular ring, reducing noise and enhancing the working environment.

Rubber bumpers inside the ring absorb impact and minimize noise levels when placing or removing the basin.

FEATURES

- Double Basin with Supporting Frame
- Supporting Frame Connected by an H-Brace
- Enhanced Strength and Stability
- Large Basin Capacity
- Unusually Quiet when Dropped into the Circular Ring

SPECIFICATIONS

- Construction: Seamless, Welded, Stainless Steel
- Stainless Steel H-Bracket
- 2" Swivel Casters
- Dimensions: 29"L x 15"W x 33"H
- Approx. Shipping Weight: 29 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)