

# Blickman 13 Qt. Lenox Kick Bucket | 7766SS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 7766SS - 13 Quart Lenox Kick Bucket has a continuous encircling rubber bumper with stainless steel cross braces that provide excellent protection for the receptacle. Constructed from heavy gauge stainless steel, this 13 qt.

pail fits perfectly into a tip-resistant base protected by rubber tips for smooth handling and noise reduction. The unit travels effortlessly on 2" swivel casters.

Limited Lifetime Warranty. Note: Click the following link to view additional models, sizes, and features available, Kick Buckets.

## FEATURES

- Continuous Encircling Rubber Bumper
- Stainless Steel Cross Braces that Provide Excellent Protection for the Receptacle
- Constructed from Heavy Gauge Stainless Steel
- Fits Perfectly into a Tip-Resistant Base
- Tip-Resistant Base Protected by rubber Tips

## SPECIFICATIONS

- Construction: Heavy Gauge Stainless Steel
- Model Number: 7766SS
- Part Number: 0817766000
- Dimensions: 16"D x 12-1/2"H
- Depth: 16"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)