

Blickman Pedestal Mounted Single-Bay Infrared Sink | 9881SSP



Stock Photo - Options May Vary

FACTORY NEW

Blickman 9881SSP - Pedestal Mounted Single-Bay Infrared Sink. The 9800 Series Lodi Scrub Sinks are built with seamless construction, with 1" covered corners. The bowl is 9-1/4" deep at the rim and 18-1/2" from spout to drain allowing for a deep scrub area.

The spout is positioned toward the rear of the sink to divert water splash away from the operator. Blickman's Pedestal Mounted Single-Bay Infrared Sink features sound-deadened material on the bowl undersides, enclosed by double-walled end panels. The sink includes an infrared sensor and knee panel for manual operation.

Sink wall mounting brackets come included (refer to the installation manual).

FEATURES

- Infra-Red Activated Single-Bay Scrub Sink
- Pedestal Mounted (wall brackets are included)
- Seamless construction
- Covered Corners
- The spout is positioned toward the rear of the sink to divert water splash away from the operator

SPECIFICATIONS

- Construction: Stainless Steel, all welding is accomplished by Heli-arc method
- Seams; seamless construction, with 1" covered corners
- Bowl Dimensions: 9-1/4" deep at the rim and 18-1/2" from spout to drain
- Infra-red sensor housing: Individual Sink
- Base: Extends 21" from the basin

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)