

# Steris Amsco Flexmatic Heavy-Duty Single Bay Scrub Sink – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Steris Amsco Flexmatic Single Bay Scrub Sink: Refurbished.

The Steris Amsco Flexmatic surgical scrub station/sink is the right answer to your pre-operative surgical scrubbing of hands and arms.

Moreover, the Steris/Amsco Flexmatic Single Bay Scrub Station/Sink Features: Knee Operation for water actuation, Stainless Steel Construction, and the formed bottom of the scrub compartment reduces splashing and flushes contaminants away from the user.

## FEATURES

- Stainless Steel construction
- Formed bottom of scrub compartment reduces splashing and flushes contaminants away from the user
- New Copper Pipes
- Upgraded Mixing Valve Kits
- New spray heads

## SPECIFICATIONS

- Overall Dimensions: 27.5"D x 32"W x 51.5"H
- Elevation from Finished Floor: 12"
- Cold Water Input: 1/2" NPT, 25-50 PSIG
- Hot Water: 1/2" NPT, 120°-140°, 25-50 PSIG
- Drain Size: 1-1/2" Overall Diameter

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)