

Blickman Wall Mounted Single-Bay Infrared Sink | 9881SSW



Stock Photo - Options May Vary

FACTORY NEW

Blickman 9881SSW - Wall Mounted Single-Bay Infrared Sink. The 9800 Series Lodi Scrub Sinks feature seamless construction with 1" covered corners. The bowl measures 9-1/4" deep at the rim and 18-1/2" from spout to drain, providing a deep scrub area.

The spout is positioned toward the rear to divert water splash away from the operator. The Blickman 9881SSW Scrub Sink features a sound-deadened underside and double-walled end panels for added durability. The sink houses a stainless steel base that will house an infra-red sensor and knee panels for manual operation.

Additionally, the housing extends 21" from the basin and hangs 14-1/2" from the finished floor. Furthermore, the Sink is supplied with wall mounting brackets.

FEATURES

- Single-Bay Sink
- Seamless construction
- Corners: Covered corners
- Activation: Infra-Red Water and Manual Knee Water/Soap Activated
- Bowl is 9-1/4" deep at rim and 18-1/2" from spout to drain allowing for a deep scrub area

SPECIFICATIONS

- Activation: Infra-Red Sensor and Knee Activated for Manual Operation
- Construction: Stainless Steel, type 300 series, nickel-bearing
- Welding: All welding is accomplished by a Heliarc method and ground smooth
- Gauges: 14 gauge
- Corner Construction: 1" covered corners

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)