

Blickman Pedestal Mounted Double-Bay Infrared Sink | 9882SSP



Stock Photo - Options May Vary

FACTORY NEW

Blickman 9882SSP - Pedestal Mounted Double-Bay Infrared Sink. The 9800 Series Lodi Scrub Sinks are built with seamless construction, with 1" covered corners. Bowls are 9-1/4" deep at the rim and 18-1/2" from spout to drain allowing for a deep scrub area.

The spout is positioned toward the rear of the sink to divert water splash away from the operator. The Blickman 9882SSP Pedestal Mounted Double-Bay Infra-Red Sink Features: Sound-deadened undersides and double-walled end panels. Each sink includes an infrared sensor and knee panels for manual operation.

Additional Features: The sink also comes with wall mounting brackets (refer to the installation manual). The sinks provide scrub areas of approximately 30"W x 21-1/2"D.

FEATURES

- Infra-Red Activated Double-Bay Scrub Sink
- Pedestal Mounted (wall brackets are included)
- Seamless construction
- Coved Corners
- The spout is positioned toward the rear of the sink to divert water splash away from the operator

SPECIFICATIONS

- Construction: Stainless Steel, all welding is accomplished by the Heli-arc method
- Seams; seamless construction, with 1" covered corners
- Bowl Dimensions (each): are 9-1/4" deep at the rim and 18-1/2" from spout to drain
- Infra-red sensor housing: Each Sink
- Base: Extends 21" from the basin

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)