

Steris Amsco Flexmatic Heavy-Duty Double Bay Scrub Sink – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Steris Amsco Flexmatic Double Bay Scrub Sink: Refurbished. The Steris Amsco Flexmatic surgical scrub station/sink is the right answer to your pre-operative surgical scrubbing of hands and arms.

Moreover, the Steris Amsco Flexmatic Double Bay Scrub Sink Features: Knee Operation for water Actuation, Stainless Steel Construction, and the formed bottom of the scrub compartment reduces splashing and flushes contaminants away from the user. Finally, they fully test each sink for leaks and proper operation before shipping.

Additionally.

FEATURES

- Stainless Steel construction
- Formed bottom of scrub compartment reduces splashing and flushes contaminants away from the user
- Upgraded Mixing Valve Kits
- New spray heads
- New On/Off Valve Diaphragm

SPECIFICATIONS

- Overall Dimensions: 27.5"D x 64"W x 51.5"H
- Elevation from Finished Floor: 12"
- Cold Water Input: 1/2" NPT, 25-50 PSIG
- Hot Water: 1/2" NPT, 120°-140°, 25-50 PSIG
- Drain Size: 1-1/2" Overall Diameter

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)