

Blickman Wall Mounted Double-Bay Infrared Sink | 9882SSW



Stock Photo - Options May Vary

FACTORY NEW

Blickman 9882SSW - Wall Mounted Double-Bay Infrared Sink. The 9800 Series Lodi Scrub Sinks feature seamless construction with 1" covered corners. The bowls measure 9-1/4" deep at the rim and 18-1/2" from spout to drain, providing a deep scrub area.

The spout positions toward the rear of the sink to divert water splash away from the operator. The Blickman 9882SSW Lodi Scrub Sink Features: The double bowls have sound-deadened undersides and double-walled end panels. Each sink includes a stainless steel base for the infrared sensor and knee panels.

Additionally, the housing extends 21" from the basin and hangs 14-1/2" from the floor. Additionally, wall mounting brackets are supplied.

FEATURES

- Double-Bay Sink
- Seamless construction
- Corners: Covered corners
- Activation: Infra-Red Water and Manual Knee Water/Soap Activated
- Bowls are 9-1/4" deep at rim and 18-1/2" from spout to drain allowing for a deep scrub area

SPECIFICATIONS

- Activation: Infra-Red Sensors and Knee Activated for Manual Operation
- Construction: Stainless Steel, type 300 series, nickel-bearing
- Welding: All welding is accomplished by a Heli-arc method and ground smooth
- Gauges: 14 gauge
- Corner Construction: 1" covered corners

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)