

# Blickman Pedestal Mounted Triple-Bay Infrared Sink | 9883SSP



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 9883SSP - Pedestal Mounted Triple-Bay Infrared Sink. Blickman's Pedestal Mounted 9800 Series Lodi Scrub Sinks feature seamless construction with 1" covered corners. The bowls measure 9-1/4" deep at the rim and 18-1/2" from spout to drain, providing a deep scrub area.

The spout positions toward the rear to divert water splash away from the operator. The Blickman 9883SSP Features: The undersides of the bowls are sound deadened and enclosed by double-walled end panels. Each sink houses an infra-red sensor and knee panels for manual operation.

Sink wall mounting brackets are included (refer to installation manual). Sinks provide scrub areas approximately 30"W x 21-1/2"D.

## FEATURES

- Infra-Red Activated Triple-Bay Scrub Sink
- Pedestal Mounted (wall brackets are included)
- Seamless construction
- Covered Corners
- The spout is positioned toward the rear of the sink to divert water splash away from the operator

## SPECIFICATIONS

- Construction: Stainless Steel, all welding is accomplished by Heli-arc method
- Seams; seamless construction, with 1" covered corners
- Bowl Dimensions (each): are 9-1/4" deep at the rim and 18-1/2" from spout to drain
- Infra-red sensor housing: Each Sink
- Base: Extends 21" from the basin

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)