

Steris Amsco Flexmatic Triple Bay Scrub Sink – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Steris Amsco Flexmatic Triple Bay Scrub Sink is a surgical scrub station designed for pre-operative hand and arm washing. It features stainless steel construction with a formed bottom compartment that reduces splashing and flushes contaminants away from the user. Knee operation provides hands-free water actuation.

This unit has been professionally refurbished to excellent operational and aesthetic condition. Upgrades include new mixing valve kits, spray heads, on/off valve diaphragms, paint, and decals. The stainless steel has been re-brushed to restore its appearance.

The sink is fully tested for leaks and operation prior to shipping. It includes a single mixing valve for each station, goose neck faucet, Z-rail, and wall brackets. A 6-month warranty is provided.

FEATURES

- Stainless Steel construction
- Formed bottom of scrub compartment reduces splashing and flushes contaminants away from the user
- Upgraded Mixing Valve Kits
- New spray heads
- New On/Off Valve Diaphragm
- Single Mixing Valve for each station
- Goose Neck Faucet
- Fully tested for leaks and operation prior to shipping

SPECIFICATIONS

- Overall Dimensions: 27.5"D x 96"W x 51.5"H
- Elevation from Finished Floor: 12"
- Cold Water Input: 1/2" NPT, 25-50 PSIG
- Hot Water: 1/2" NPT, 120°-140°, 25-50 PSIG
- Drain Size: 1-1/2" Overall Diameter
- Electricity (if using optional infrared): 120V, 1 Amp

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)