

Blickman Wall Mounted Triple-Bay Infrared Sink | 9883SSW



Stock Photo - Options May Vary

FACTORY NEW

Blickman 9883SSW - Wall Mounted Triple-Bay Infrared Sink. The 9800 Series Lodi Scrub Sinks feature seamless construction with 1" covered corners. The bowls measure 9-1/4" deep at the rim and 18-1/2" from spout to drain, providing a deep scrub area.

The design positions the spout toward the rear of the sink to divert water splash away from the operator. The Blickman 9883SSW Lodi Scrub Sink Features: Undersides of the triple bowls are sound deadened and enclosed by double-walled end panels. Each sink houses a stainless steel base that will house an infra-red sensor and knee panels for manual operation.

The housing extends 21" from the basin and hangs 14-1/2" from the finished floor.

FEATURES

- Triple-Bay Sink
- Seamless construction
- Corners: Coved corners
- Activation: Infra-Red Water and Manual Knee Water/Soap Activated
- Bowls are 9-1/4" deep at rim and 18-1/2" from spout to drain allowing for a deep scrub area

SPECIFICATIONS

- Activation: Infra-Red Sensors and Knee Activated for Manual Operation
- Construction: Stainless Steel, type 300 series, nickel-bearing
- Welding: All welding is accomplished by a Heli-arc method and ground smooth
- Gauges: 14 gauge
- Corner Construction: 1" coved corners

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)